

OBESITY DETECTION USING MACHINE LEARNING

Step-by-Step Flow

1. Anaconda Installation steps

Step 1- Download Anaconda: Visit the official Anaconda website

Install the latest version of Anaconda

<https://www.filehorse.com/download-anaconda/>

Next-Next-All Users-Install

Step 2-

1. **Install Anaconda:** Run the installer and follow the on-screen instructions. Choose the recommended settings unless you have specific preferences.
2. **Create a Virtual Environment:** Once Anaconda is installed, open the Anaconda Navigator or Anaconda Prompt.
 - Create a virtual environment for your project using the command:
conda create -n Obesity Detection python=3.8
 - Create environment 3.8.17

Step 3- Select environment 3.8.17 then Spyder install and launch



Rename Envo. Obesity Detection



(Right Click) Open terminal



Install Lib

- **Libraries**

1] pip install tensorflow

2] pip install matplotlib

3] pip install keras

4] pip install pillow

5] pip install numpy

6] pip install opencv-python

7] pip install matplotlib

8] pip install scikit-learn

9] pip install tkvideo

10] pip install mediapipe

11] pip install gts

12] pip install pandas

13] pip install flask

14] pip install mixtend

2.Install DB Browser (SQLite)

Step 1: Download SQLite: Visit the official SQLite website

Install the latest version of SQLite

<https://sqlitebrowser.org/dl/>

Step 2: Connect to the MySQL Database

- Ensure your MySQL server is running.
- Open SQLite -> New Database -> Open Folder (OBD-2024 CODE) -> Select File (Evaluation, registration).
- Go To Execute SQL, And Use Command: **(SELECT * FROM REGISTRATION)**, To Retrieve Data.

3.Spyder

Step 1:

In Anaconda Navigator, launch Spyder within the Obesity detection environment.

Step 2:

- Open the Project Folder in Spyder
- In Spyder, open r project folder using: **File -> Open (Folder) ->"OBD-2024 CODE"**.
- Run File **"Gui main.py"** and Run It.